

BRENNEN SKIN GRAFT MESHER

*The Non-Crushing, Non-Carrier Based Instrument
That Makes Meshing Safe and Easy.*

The Brennen Skin Graft Mesher is a totally self-contained, state-of-the-art meshing mechanism. The instrument meshes by "pinching" not cutting the graft, creating a uniform pattern and maintaining graft integrity.



Carrier-Free

The Brennen Skin Graft Mesher does not require plastic carriers, allowing you to mesh any length graft without being dictated by the length of a carrier. This also saves the time, trouble and cost of ordering disposable carriers.

Graft Versatility

To meet your meshing needs, the Brennen Skin Graft Mesher is available in two platform widths (4.5" and 6.5") and a variety of expansion ratios. Brennen Skin Graft Meshers accept any length graft and cope well with both thin and thick grafts¹.

Safe and Easy to Use

The Brennen Skin Graft Mesher is designed for convenience and ease of use. The smooth operation of the ratchet handle makes the meshing process effortless. The weight of the base provides added stability and holds the mesher firmly in place for unassisted operation.

The Brennen Skin Graft Mesher does not contain blades, eliminating the need for sharpening and reducing the concern for injury. The crossbar gives an extra element of steadiness when meshing and makes the unit easily transportable.



BRENNEN
MEDICAL

www.brennenmed.com

reliable performance. | naturally
extraordinary outcomes. |

BRENNEN SKIN GRAFT MESHER

Ordering Information

4.5 inch Brennen Skin Graft Mesher

meshes grafts up to 4.5 in (11.5 cm) wide

CATALOG NO.	DESCRIPTION (APPROX.)
BSGM-PC	Pie Crust (minimal perforations)
BSGM-0011	1:1 Expansion Ratio
BSGM-0021	2:1 Expansion Ratio
BSGM-0031	3:1 Expansion Ratio
BSGM-0041	4:1 Expansion Ratio
BSGM-0061	6:1 Expansion Ratio
BSGM-0100	Mesher Storage/Sterilizer Box

6.5 inch Brennen Skin Graft Mesher

meshes grafts up to 6.5 in (16.5 cm) wide

CATALOG NO.	DESCRIPTION (APPROX.)
BSGX-0011	1:1 Expansion Ratio
BSGX-0021	2:1 Expansion Ratio
BSGX-0031	3:1 Expansion Ratio
BSGX-0041	4:1 Expansion Ratio
BSGX-0061	6:1 Expansion Ratio

Satisfaction Guarantee

Your complete satisfaction with the Brennen Skin Graft Mesher is guaranteed. Brennen Medical will refund the full purchase price within 30 days of purchase date if you are not totally convinced the Brennen Skin Graft Mesher meets your needs.*

A free 30-day evaluation of a Brennen Skin Graft Mesher is available.* Contact your local sales representative or call customer service at (800) 328-9105 for more information.

**For More Information or
To Order Call (800) 328-9105**

Brennen Medical, LLC
1290 Hammond Road
Saint Paul, MN 55110
Ph: (651) 429-7413 • (800) 328-9105
Fax: (651) 429-8020
100616-10/09

Pie Crust
(minimal perforations)



1:1



2:1



3:1



4:1



6:1



Representation of Brennen Skin Graft Mesher patterns using paper model



BRENNEN
MEDICAL

¹ Taghizadeh R, Gilbert PM. Comparison of commonly used mesher types in burn surgery revisited. Burns 2008;34:109-110.

* Available in the U.S. only