



# Bone Instruments



# Bone Instruments

**Brauer Rocking Nasal Chisel**

Catalog#	Length	Width
P708	6 1/4" (159 mm)	13 mm



P708

**Sheen Acute Angle Radix Chisel**

Catalog#	Length	Width
P883	7" (178 mm)	10 mm



P883

**McIndoe Nasal Chisel**

Catalog#	Length	Width
P15869	5 1/2" (140 mm)	9 mm
P158611	5 1/2" (140 mm)	11 mm
P158613	5 1/2" (140 mm)	13 mm
P158615	5 1/2" (140 mm)	15 mm




P158613

**McIndoe Nasal Chisel with Mushroom Head**

Catalog#	Length	Width
P15929	5 1/2" (140 mm)	9 mm
P159211	5 1/2" (140 mm)	11 mm
P159213	5 1/2" (140 mm)	13 mm
P159215	5 1/2" (140 mm)	15 mm



P159213



## Braithwaite's Nasal Chisel

Catalog#	Description	Length	Width
P1554A	Right	5 1/2" (140 mm)	5 mm
P1554B	Left	5 1/2" (140 mm)	5 mm



## Silver Cartilage Chisel

P1554A


Catalog#	Description	Length	Width
P1630A	Straight	7 1/4" (184 mm)	7 mm
P1630B	Curved Left	7 1/4" (184 mm)	7 mm
P1630C	Curved Right	7 1/4" (184 mm)	7 mm



## Freer Septum Chisel

P1630C

Catalog#	Description	Length	Width
P4414	Straight	6 1/2" (165 mm)	4 mm
P4415	Curved	6 1/2" (165 mm)	4 mm
P4412	Straight	7 1/2" (191 mm)	4 mm
P4413	Curved	7 1/2" (191 mm)	4 mm



## Clodius Bone Chisel

P4414

Catalog#	Length	Width
P3813	8 3/4" (222 mm)	6 mm
P3814	8 3/4" (222 mm)	8 mm
P3815	8 3/4" (222 mm)	10 mm
P3817	8 3/4" (222 mm)	16 mm
P3818	8 3/4" (222 mm)	20 mm



P3814

**Fomon Chisel**

Catalog#	Description	Length	Width
P7500	Guarded Point	6 1/2" (165 mm)	5 mm
P7502	Guarded Point	6 1/2" (165 mm)	7 mm
P7503	Guarded Point	6 1/2" (165 mm)	8 mm



P7502

**Cottle Chisel**

Catalog#	Length	Width
P7510	7 1/4" (184 mm)	4 mm
P7515	7 1/4" (184 mm)	7 mm
P7517	7 1/4" (184 mm)	9 mm
P7520	7 1/4" (184 mm)	12 mm



P7510

**Cottle Chisel**

Catalog#	Description	Length	Width
P7530	Curved	7 1/4" (184 mm)	6 mm



P7530

**Hajek Septum Chisel**

Catalog#	Description	Length	Width
P7570	Small	6 1/4" (159 mm)	4 mm
P7572	Medium	6 1/4" (159 mm)	6 mm
P7574	Large	6 1/4" (159 mm)	8 mm



P7570



## Sheehan Chisel

Catalog#	Length	Width
P7600	6" (152 mm)	7 mm
P7602	6" (152 mm)	8 mm
P7604	6" (152 mm)	10 mm
P7606	6" (152 mm)	11 mm



P7600



## Rees Nasal Osteotome\*

Catalog#	Length	Width
P711	7 3/4" (197 mm)	3 mm
P712	7 3/4" (197 mm)	4 mm

\*Anodized Aluminum handle with stainless steel working end.



P711

## French Pattern Osteotome

Catalog#	Length	Width
P15025	5 3/4" (146 mm)	5 mm
P15027	5 3/4" (146 mm)	7 mm
P15029	5 3/4" (146 mm)	9 mm




P15025

## Rozner Nasal Osteotome

Catalog#	Length	Width
P1520	6 1/2" (165 mm)	6.5 mm



P1520



### Ward Nasal Osteotome

Catalog#	Length	Width
P15463	6 1/2" (165 mm)	3 mm
P15465	6 1/2" (165 mm)	5 mm



P15463

### Lambotte Osteotome

Catalog#	Length	Width
P3800	5" (127 mm)	4 mm
P3802	5" (127 mm)	6 mm
P3804	5" (127 mm)	8 mm
P3806	5" (127 mm)	10 mm
P3808	5" (127 mm)	12 mm



P3800

### Neivert Osteotome

Catalog#	Description	Length	Width
P4761	Single Guarded, Straight	8" (203 mm)	5 mm
P4885	Curved, Single Guarded Right	8" (203 mm)	4 mm
P4886	Curved, Single Guarded Left	8" (203 mm)	4 mm



P4761

**Cinelli Osteotome**

Catalog#	Description	Length	Width
P4762	Single Guarded	6 1/4" (159 mm)	13 mm


**Cinelli Osteotome**

P4762

Catalog#	Description	Length	Width
P4763	Double Guarded	6 1/4" (159 mm)	10 mm
P4786	Double Guarded	6 1/4" (159 mm)	12 mm
P4764	Double Guarded	6 1/4" (159 mm)	14 mm
P4765	Double Guarded	6 1/4" (159 mm)	16 mm
P4983	Double Guarded, Curved	6 1/4" (159 mm)	10 mm



**Padgett Osteotome**

P4763

Catalog#	Length	Width
P4810	7" (178 mm)	2 mm
P4811	7" (178 mm)	3 mm
P4812	7" (178 mm)	4 mm
P4813	7" (178 mm)	5 mm
P4814	7" (178 mm)	6 mm
P4815	7" (178 mm)	7 mm
P4816	7" (178 mm)	8 mm
P4817	7" (178 mm)	9 mm
P4818	7" (178 mm)	10 mm
P4819	7" (178 mm)	11 mm
P4820	7" (178 mm)	12 mm
P4828	7" (178 mm)	15 mm
P4830	Set of 12 Padgett Osteotomes	



P4810



### Tessier Maxillofacial Osteotome

Catalog#	Description	Length	Width
P4935	Straight	8 1/2" (216 mm)	10 mm
P4936	Straight	8 1/2" (216 mm)	15 mm
P4940	Curved	8 1/2" (216 mm)	10 mm
P4941	Curved	8 1/2" (216 mm)	15 mm



### Planas T Osteotome

Catalog#	Description	Length	Width
P4979	Double-Guarded, Curved	7 1/2" (191 mm)	5 mm
P4985	Double-Guarded, Curved	7 1/2" (191 mm)	8 mm

P4941



P4979



P4985

### Giunta Double-Guarded Osteotome

Catalog#	Description	Length	Width
P1632	Curved	6 1/4" (159 mm)	5 mm
P1634	Straight	6 1/4" (159 mm)	5 mm



### Furnas Bayonet Osteotome

Catalog#	Length	Width
P4981	8" (203 mm)	13 mm

P1632



P4981

**Rish Osteotome**

Catalog#	Description	Length	Width
P4982	Guarded	7" (178 mm)	6.5 mm



P4982

**Rubin Nasal Osteotome**

Catalog#	Description	Length	Width
P7525	With Stabilizer	6 1/2" (165 mm)	11 mm
P7526	With Stabilizer	6 1/2" (165 mm)	14 mm



P7525

**Pitanguy Nasal Osteotome**

Catalog#	Length	Width
P7535	8" (203 mm)	8 mm
P7536	8" (203 mm)	14 mm
P7537	8" (203 mm)	16 mm



P7535



## Freer Gouge


Catalog#	Description	Length	Width
P4417	Straight	6 1/4" (159 mm)	4 mm
P4418	Straight	6 1/4" (159 mm)	6 mm



P4417

## Heden Bone Gouge

Catalog#	Description	Length	Tip Width
P7670	Straight	8 1/2" (216 mm)	10 mm



P7670

## Killian Septum Gouge

Catalog#	Length	Width
P4419	6 1/4" (159 mm)	5 mm



P4419

**West Bone Gouge**

Catalog#	Length	Width
P7665	7 1/2" (191 mm)	6 mm




P7665

**Parkes Hump Gouge**

Catalog#	Description	Length	Width
P7595	Double Guarded	7" (178 mm)	6 mm



P7595



### Blumenthal Bone Rongeur

Catalog#	Description	Length	Bite
P4462	Curved	6" (152 mm)	4 mm



P4462

### Blumenthal Bone Rongeur

Catalog#	Description	Length	Bite
P4463	Narrow	7 1/4" (184 mm)	3 mm




P4463

### Beyer Rongeur

Catalog#	Description	Length	Bite
P4465	Delicate	7" (178 mm)	3 mm



P4465



## Hardee Bone Rongeur

Catalog#	Description	Length	Bite
P4868	Slight Curve	6" (152 mm)	2 mm
P4869	Curved	6" (152 mm)	2 mm
P4784	Curved	6" (152 mm)	4 mm



P4869

## Kleinert-Kutz Rongeur


Catalog#	Description	Length	Bite
P8140	Strong Curve	5 3/4" (146 mm)	3 mm




P8140

## Friedman Rongeur

Catalog#	Description	Length	Bite
P4775	Delicate Jaw	5" (127 mm)	3 mm



P4775



**Prentis Rongeur**

Catalog#	Length	Bite
P4776	5 1/2" (140 mm)	2.5 mm



P4776

**Ruskin Bone Rongeur**

Catalog#	Description	Length	Bite
P8175	Curved	7 1/4" (184 mm)	5 mm




P8175

**Ferris-Smith Spurling Rongeur**

Catalog#	Length	Bite
P4468	7" (178 mm)	4 mm x 10 mm



P4468



## Cottle-Kazanjan Nasal Hump Forceps

Catalog#	Description	Length
P8125	Narrow Jaws, Straight	7 1/2" (191 mm)



P8125

## Rowland Nasal Hump Forceps

Catalog#	Description	Length
P8130	Double Action, Straight Jaw	7" (178 mm)




P8130

## Liston Bone Cutting Forceps

Catalog#	Description	Length
P1418	Straight	7 1/2" (191 mm)
P1420	Angled-on-Flat	7 1/2" (191 mm)



P1418



## McIndoe Bone Cutting Forceps

Catalog#	Description	Length
P1422	Angled-on-Flat, Double Action Jaws	7 1/4" (184 mm)



P1422

## Padgett Bone Cutting Forceps

Catalog#	Description	Length
P4268	Straight	5 3/4" (146 mm)
P4269	Angled	5 3/4" (146 mm)



P4268

## Key Bone Cutting Forceps

Catalog#	Description	Length
P1414	Double Action Joints	10" (254 mm)



P1414

## Ruskin-Liston Bone Forceps

Catalog#	Description	Length
P4458	Wide, Straight	7 1/4" (184 mm)
P4459	Wide, Angled	7 1/4" (184 mm)



P4459

**Volkman Scoop**

Catalog#	Description	Length	Scoop 1	Scoop 2
P1476A	Double-Ended	8 3/4" (222 mm)	3 mm x 5 mm	5 mm x 10 mm
P1476B	Double-Ended	8 3/4" (222 mm)	4 mm x 7 mm	6 mm x 12 mm
P1476C	Double-Ended	8 3/4" (222 mm)	5 mm x 9 mm	7 mm x 13 mm




P1476A

**Meyhoeffer Curette**

Catalog#	Cup Size	Length	Cup Diameter
P7989	Size 000	5" (127 mm)	0.5 mm
P7990	Size 00	5" (127 mm)	1 mm
P7991	Size 0	5" (127 mm)	1.5 mm
P7992	Size 1	5" (127 mm)	1.75 mm
P7993	Size 2	5" (127 mm)	2 mm
P7994	Size 3	5" (127 mm)	2.5 mm
P7995	Size 4	5" (127 mm)	3.5 mm
P7996	Size 5	5" (127 mm)	4 mm
P7997	Size 6	5" (127 mm)	4.5 mm



P7995



**Fox Dermal Curette**

Catalog#	Length	Tip Diameter
P7975	5 1/2" (142 mm)	3 mm
P7976	5 1/2" (142 mm)	4 mm
P7977	5 1/2" (142 mm)	5 mm
P7978	5 1/2" (142 mm)	6 mm



P7975

