



Mallets and Cartilage Instruments



Mallets and Cartilage Instruments

Fowler Aluminum Mallet

Catalog#	Description	Length	Diameter
P805	10 oz.	7 3/4" (197 mm)	1 1/2" (38 mm)



P805

Lucae Mallet

Catalog#	Description	Length	Diameter
P4215	Stainless Steel	8" (203 mm)	1" (25 mm)



P4215

Mallet

Catalog#	Description	Length	Diameter
P1492	Solid Stainless Steel, 3/4 lb.	7 3/4" (197 mm)	35 mm



P1492



Quisling Intranasal Hammer

Catalog# Length

P3780 7" (178 mm)



P3780

Gerzog Mallet

Catalog# Description Length Diameter

P4389 Chrome Plated 6 3/4" (171 mm) 24 mm



P4389

Kirk Mallet

Catalog# Length Diameter

P3779 7 3/4" (197 mm) 41 mm




P3779

Padgett Nasal and Orthopedic Mallet


Catalog# Description Length Diameter

P4398 With Replaceable Bakelite Faces 7 1/2" (191 mm) 1" (25 mm)

P4398R Replacement Bakelite Faces/Pair



P4398



Mead Mallet

Catalog#	Description	Length	Diameter
P4399	Nylon Faces	7 1/2" (191 mm)	20 mm



P4399

Padgett Nasal Mallet

Catalog#	Description	Length	Diameter
P4984	Nylon Faces	7 3/8" (187 mm)	1" (25 mm)
P4984R	Insert Replacement For P4984/Pair		1" (25 mm)



P4984



P4984R

Sherman Nasal Mallet

Catalog#	Description	Length	Diameter
P4900	Copper Head	7 1/8" (182 mm)	3/4" (19 mm)



P4900



260
250
240
230
220
210
200
190
180
170
160
150
140
130
120
110
100
90
80
70
60
50
40
30
20
10

Cottle Cartilage Crusher


Catalog#	Description	Length	Width
P666	Smooth Bed	3" (76 mm)	16 mm
P667	Serrated Bed	3" (76 mm)	16 mm



P666

Jost Cartilage Crusher


Catalog#	Bed
P3575	64 mm x 32 mm



P3575

Lower Lateral Cartilage Morselizer

Catalog#	Description	Length	Jaw
P3819	Serrated Jaws with Pyramid Shaped Teeth	7 1/2" (191 mm)	4.5 mm x 9.5 mm



Rubin Morselizer

Catalog#	Description	Length	Jaw
P3820	Rubin Morselizer with Clip	7 1/2" (191 mm)	11 mm x 6 mm
P3829	Clip for P3820		




Padgett-Rubin Morselizer

Catalog#	Description	Length	Jaw
P3830	Angled Handles, Serrated Long Jaws with Pyramid Shaped Teeth and Clip	8 1/4" (210 mm)	15 mm x 5 mm
P3831	Replacement Clip for P3830		



P3830



Gorney Septal Morselizer


Catalog#	Description	Length
P2215	Double Action, Angled, Tungsten Carbide	8 1/4" (210 mm)



P2215

Gunter Lower Lateral Morselizer

Catalog#	Description	Length
P2216	Tungsten Carbide Serrated Jaws, 1-Smooth Jaw and 1-Serrated Jaw with Pyramid Shaped Teeth	7" (178 mm)
P2211	Tungsten Carbide Jaws with Semi-Sharp Longitudinal Ridges	7" (178 mm)



P2216