




Skin Graft Products




Skin Graft Products




Slim Design, Serious Performance
Model S Slimline Electric Dermatome Kit

No.	Catalog#	Description
1.	3539700	Model S Slimline Dermatome Set, Complete
2.	DP0009	Model S Slimline Dermatome Handpiece ☆ ▲
3.	3535399100	Head Assembly Guard Plate ☆
4.	35392661 35392662 35392663 35392664	1" (2.54 cm) Width Blade Clip (Not Shown) ▲ 2" (5.08 cm) Width Blade Clip ☆ 3" (7.62 cm) Width Blade Clip ☆ 4" (10.16 cm) Width Blade Clip ☆
5.	353971513 353971525	13' (3.96 m) Handpiece Power Cord ☆ 25' (7.62 m) Handpiece Power Cord (Not Shown) ▲
6.	DP0007	Universal Power Supply ☆ 100–250 Volts
7.	3539620	U.S. Power Supply Wall Cord ☆ (Specify Country of Use for Alternative Plugs)
8.	3539259	Calibration Guide ☆
9.	3539257	Screwdriver for Dermatome, Philips, Stainless Steel ☆
10.	3539240	Carrying Case, ABS plastic (Not Shown) ☆
11.	3539252	Sterile Dermatome Blades, Box of 10 (Not included in Kit)

☆ Product Is Included in Its Complete Kit
▲ Special Order, if ordered separately






Model S6 Slimline Electric Dermatome Kit

No.	Catalog#	Description
1.	3539900	Model S6 Slimline Dermatome Set, Complete
2.	DP0010	Model S6 Slimline Dermatome Handpiece ☆ ▲
3.	35399100	Head Assembly Plate (Not Shown) ☆
4.	35399101	1" (2.54 cm) Width Blade Clip (Not Shown) ▲
	35399102	2" (5.08 cm) Width Blade Clip ☆
	35399103	3" (7.62 cm) Width Blade Clip ☆
	35399104	4" (10.16 cm) Width Blade Clip ☆
	35399105	5" (12.7 cm) Width Blade Clip ☆
	35399106	6" (15.2 cm) Width Blade Clip ☆
5.	353971513	13' (3.96 m) Handpiece Power Cord ☆
	353971525	25' (7.62 m) Handpiece Power Cord (Not Shown) ▲
6.	DP0007	Universal Power Supply ☆ 100–250 Volts
7.	3539620	U.S. Power Supply Wall Cord ☆ (Specify Country of Use for Alternative Plugs)
8.	3539259	Calibration Guide ☆
9.	3539257	Screwdriver for Dermatome, Philips, Stainless Steel ☆
10.	3539240	Carrying Case, ABS plastic (Not Shown) ☆
11.	3539602	Sterile Dermatome Blades, Box of 10 (Not Shown) (Not included in Kit)


☆ Product Is Included in Its Complete Kit
▲ Special Order, if ordered separately




Model SB Slimline Electric Dermatome Kit

No.	Catalog#	Description
1.	3539800	Model SB Slimline Dermatome Set, Complete
2.	DP0810	Model SB Slimline Dermatome Handpiece ☆ ▲
3.	35392660	Head Assembly Plate (Not Shown) ☆
4.	35392661 35392662 35392663 35392664	1" (2.54 cm) Width Blade Clip (Not Shown) ▲ 2" (5.08 cm) Width Blade Clip ☆ 3" (7.62 cm) Width Blade Clip ☆ 4" (10.16 cm) Width Blade Clip ☆
5.	353971513 353971525	13' (3.96 m) Handpiece Power Cord ☆ 25' (7.62 m) Handpiece Power Cord (Not Shown) ▲
6.	DP0800	Universal Power-Supply 100-250 volts ☆
7.	3539620	U.S. Power Supply Wall Cord ☆ (Specify Country of Use for Alternative Plugs)
8.	3539259	Calibration Guide ☆
9.	3539257	Screwdriver for Dermatome, Philips, Stainless Steel ☆
10.	3539240	Carrying Case, ABS plastic (Not Shown) ☆
11.	3539252	Sterile Dermatome Blades, Box of 10 (Not Shown) (Not included in Kit)

☆ Product Is Included in Its Complete Kit
▲ Special Order, if ordered separately






Model PI Electric Dermatome Kit

No.	Catalog#	Description
1.	3539600 3539600230	Model PI Dermatome Set, Complete, 120 Volts Model PI Dermatome Set, Complete, 230 Volts (Not Shown)
2.	N/A	Model PI Dermatome Handpiece and Cord ☆ ▲
3.	35396100	Head Assembly Guard Plate ☆
4.	35396101 35396102 35396103 35396104 35396105 35396106	1" (2.54 cm) Width Blade Clip (Not Shown) ▲ 2" (5.08 cm) Width Blade Clip ☆ 3" (7.62 cm) Width Blade Clip ☆ 4" (10.16 cm) Width Blade Clip ☆ 5" (12.7 cm) Width Blade Clip ☆ 6" (15.24 cm) Width Blade Clip ☆
5.	DP0005 DP0006 3539620	Universal Power Supply (120 Volts) ☆ Universal Power Supply (230 Volts) (Not Shown) ☆ U.S. Power Supply Wall Cord ☆ (Specify Country of Use for Alternative Plugs)
6.	3539259	Calibration Guide ☆
7.	3539257	Screwdriver for Dermatome, Philips, Stainless Steel ☆
8.	3539240	Carrying Case, ABS Plastic (Not Shown) ☆
9.	3539602	Sterile Dermatome Blades, Box of 10 (Not Shown) (Not Included in Kit) ☆

☆ Product Is Included in Its Complete Kit

▲ Special Order, if ordered separately




Model B Electric Dermatome Kit

No.	Catalog#	Description
1.	3539250 3539250230	Model B Dermatome Set, Complete, 120 Volts Model B Dermatome Set, Complete, 230 Volts (Not Shown)
2.	DP0001	Model B Dermatome Handpiece and Cord ☆ ▲
3.	35392550	Head Assembly Guard Plate ☆
4.	35392551 35392552 35392553 35392554	1" (2.54 cm) Width Clip (Not Shown) ▲ 2" (5.08 cm) Width Clip ☆ 3" (7.62 cm) Width Clip ☆ 4" (10.16 cm) Width Clip ☆
5.	DP0002 DP0004	Universal Power Supply (120 Volts) ☆ Universal Power Supply (220 Volts) ☆ (Not Shown)
6.	3539259	Calibration Guide ☆
7.	3539257	Screwdriver for Dermatome, Philips, Stainless Steel ☆
8.	3539240	Carrying Case, ABS Plastic (Not Shown) ☆
9.	3539252	Sterile Dermatome Blades, Box of 10 (Not Shown) (Not Included in Kit) ☆

☆ Product Is Included in Its Complete Kit

▲ Special Order, if ordered separately





Model C Air Dermatome Kit

No.	Catalog#	Description
1.	3539500	Model C Air Dermatome Set, Complete
2.	AP0001	Model C Air Dermatome Handpiece ☆ ▲
3.	35392660	Head Assembly Guard Plate ☆
4.	35392661 35392662 35392663 35392664	1" (2.54 cm) Width Blade Clip (Not Shown) ▲ 2" (5.08 cm) Width Blade Clip ☆ 3" (7.62 cm) Width Blade Clip ☆ 4" (10.16 cm) Width Blade Clip ☆
5.	3539401	12" (3.66 m) Air Hose and Connector for Model C Dermatome ☆
6.	3539259	Calibration Guide ☆
7.	3539257	Screwdriver for Dermatome, Philips, Stainless Steel ☆
8.	3539240	Carrying Case, ABS Plastic (Not Shown) ☆
9.	3539252	Sterile Dermatome Blades, Box of 10 (Not Shown) (Not Included in Kit) ☆


☆ Product Is Included in Its Complete Kit

▲ Special Order, if ordered separately

Powered Dermatome Accessories

Catalog# Description

3539799 Dermatome Sterilization Tray



3539799

353971525 Model S Removable Handpiece Cord, 25' (7.62 m)



353971525

3539257 Screwdriver for Dermatomes, Phillips, Stainless Steel



3539257

3539259 Calibration Guide, Stainless Steel



3539259

3539701 Dermatome Hex Nuts


(Picture Not Shown)

3539701

3539702 Dermatome Wrench



3539702






Manual Drum Dermatome Kit

No.	Catalog#	Description
1.	3539100TA 3539100BTA	Regular Dermatome Complete with Accessories and Single Use Blades Baby Dermatome Complete with Accessories and Single Use Blades (Not Shown)
2.	3539101 3539102 3539105 3539106	Dermatome (Drum and Stand Only) 4x8 inch Drum, Regular Size ✦ Baby Size Dermatome (Drum and Stand Only) 3x8 inch (Not Shown) + Sterile, Single Use Blades, Regular Size, Box of 10 (Not Shown) (2 Boxes of Blades Included with Initial Order) ✧ Sterile, Single Use Blades, Baby Size, Box of 10 (Not Shown) (2 Boxes of Blades Included with Initial Order) ✧
	3539107*	Adapter for Use with Regular Size Single Use Blades (Not Shown)
	3539108*	Adapter for Use with Baby Size Single Use Blades (Not Shown)
	3539109	Honing Handle, Regular Size (Not Shown) ✧
	3539110	Honing Handle, Baby Size (Not Shown)
	3539111*	Adapter for Honing Handle, Regular Size, for Use with Regular Single Use Blades (Not Shown) ✧
	3539112*	Adapter for Honing Handle, Baby Size, for Use with Baby Single Use Blades (Not Shown) ✧
	3539115	Carrying Case (Leatherette) (Not Shown) ✧
3.	3539118	Clip for Replacement; for Holding Blades on Blade Bar ✧
4.	3539117 35391520	Thickness Gauge ✧ PADGETT Dermatome Tape, 10/box (Not Shown)
5.	3539102	Drum Stand (Smaller Stand for Baby Dermatome) ✧

✧ Product Is Included in Its Complete Kit
 ✦ uses No. 105 Regular Size Blades
 + uses No. 106 Baby Size Blades

* These products are no longer necessary on newer versions of the product or after a repair of the blade bar. The newer version serial numbers are N7910 and higher, and 5IN780 and higher.




Brennen Skin Graft Mesher System

Catalog#	Approx. Mesh Ratio Guide
P902701	1:1
P902702	2:1
P902703	3:1
P902704	4:1
P902706	6:1
P902708	8:1
P9027SC	Mesher Carrying Case

System includes tool kit: Allen Wrench, Screwdriver, Mineral Oil, and Soft Bristled Brush.

* This product available only in the United States.




Skin Graft Expanders

Catalog#	Description
P170	Mesh Skin Graft Expander, 4 1/4" x 5 1/2"



P170

P172	Teflon Roller for Skin Graft Expander
------	---------------------------------------



P172

P173	Roller and Handle, Complete for P170
------	--------------------------------------



P173

Catalog# Description

P8632 USF Debriding Forceps



P8632

P4597 WATSONS MOD BRAITHWAITE Skin Graft Knife, 12" Long. Uses a 6 1/4" Blade
P4597LH WATSONS MOD BRAITHWAITE Skin Graft Knife, 12" Long. Uses a 6 1/4" Blade
Left Handed




P4597

P4598 COBBETTS MOD BRAITHWAITE Skin Graft Knife, 12" Long. Uses a 6 1/4" Blade
P4596 Sterile Blades for BRAITHWAITE Skin Graft Knife, 6 1/4" Long, 10/Package for
P4597, P4597LH, P4598



P4598

P14701 SILVERS Miniature Skin Graft Knife, 7 1/2". Uses Double-Edge Razor Blade



P14701

