

S C A N L A N[®]

Fast-F.I.C.S.[™]

FORCEPS AND INSTRUMENT CLEANING STATION



SCANLAN ★[®]

INTERNATIONAL

Products Designed and Manufactured by the Scanlan Family Since 1921

S C A N L A N[®]

Fast-F.I.C.S.[™]

FORCEPS AND INSTRUMENT CLEANING STATION


The SCANLAN[®] Fast-F.I.C.S.[™] forceps and instrument cleaning station is a single-use product designed to provide fast, easy and effective removal of tissue and debris which commonly gathers on the tips of delicate forceps, haemostats, clamps and other surgical instrumentation.

This compact, 3-in-1 station offers multi-layered cleaning bristles, a non-abrasive cleaning cloth and a cautery cleaning pad — all mounted on a peel-off adhesive backing for secure placement near the surgical field.

Unique and easy to use, the SCANLAN[®] Fast-F.I.C.S.[™] forceps and instrument cleaning station reduces time, costs, and inconvenience within the operating room.

The SCANLAN[®] Fast-F.I.C.S.[™] forceps and instrument cleaning station features...

- Multi-layered cleaning bristles effectively remove tissue and debris with “one pass”
- A non-abrasive cleaning cloth for quick secondary cleaning of delicate forceps and instrument tips
- A cautery pad cleans uncoated cautery blades thus eliminating the need for additional cleaning devices
- An adhesive backing for fast and secure placement within the sterile field or on the instrument stand
- Compact design allows for maximum utilization in the sterile field
- Latex safe
- Sterile packaged
- Single-use
- Radiopaque



The SCANLAN[®] Fast-F.I.C.S.[™] forceps and instrument cleaning station was designed in conjunction with renowned vascular surgeon Russell H. Samson, M.D.

ORDERING INFORMATION

1001-001 SCANLAN[®] Fast-F.I.C.S.[™] forceps and instrument cleaning station
1 per sterile package, 5 packages per box

**For more information or to order,
contact your local Scanlan representative or call...**

U.S. & Canada...800-328-9458

International...651-298-0997

FAX...651-298-0018



INTERNATIONAL

Products Designed and Manufactured by the Scanlan Family Since 1921

ONE SCANLAN PLAZA • SAINT PAUL, MINNESOTA 55107 • USA
800-328-9458 • 651-298-0997 • FAX 651-298-0018 • www.scanlangroup.com

© 1999 Scanlan International, Inc.
MK-B-016 Rev A