

PINCAPTM

Orthopedic Pin Protectors

- Fits all common pins and wires
- Protects patient and OR professionals
- Double wall design
- Color coded
- Sterile packaged
- Radiopaque
- Latex safe



(Photo enlarged to show details)



INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921

PINCAPTM

Orthopedic Pin Protectors

A unique single-use protective cover for exposed ends of orthopedic pins and wires. PinCapTM orthopedic pin protectors will help prevent the following:

- Dermal abrasions
- Damage to bedding
- Clothing snags
- Surgical glove snags
- Annoying problems associated with exposed pins

Features

- Radiopaque
- Double-wall design
- Color coded
- Sterile packaged:
2 per package / 20 packages per box

Benefits

- Fits all common pins and wires
- Patient safety
- Added protection
- Easily identifiable
- Ease to use

Pin or wire size							
Product No.	Color	Inches			mm		
301-0110	White (fits 3 sizes K-wire)	.028	.035	.045	.71	.89	1.14
301-0120	Pink & White	.062				1.57	
301-0130	Blue & White	.078		5/64		1.98	

For more information or to place an order, please contact your local Scanlan representative, or call...

U.S. & Canada: 800-328-9458

International: 651-298-0997

FAX: 651-298-0018 • www.scanlangroup.com



INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921
ONE SCANLAN PLAZA • SAINT PAUL, MINNESOTA 55107 • U.S.A.

MK-B-Pincap Rev A
© 2006 Scanlan International, Inc.