

INSTRUMENT CARE PRODUCTS ■ Sterilization Trays ■ Metal



■ SCANLAN® Metal Sterilization Trays, Aluminum, Silicone Inserts, Detachable Lids

Cat. No.	Description	Instrument Capacity
NEW * 2138-03	Silver-gray color, 10" x 10" x 1.5" (25 cm x 25 cm x 3.8 cm)	8-10
NEW * 2140-03	Silver-gray color, 10" x 15" x 1.5" (25 cm x 38 cm x 3.8 cm)	10-12
NEW * 2141-03	Silver-gray color, base with accessory compartment, additional tray insert, 10" x 15" x 2 5/8" (25 cm x 38 cm x 6.6 cm)	15+
NEW * 2142-03	Silver-gray color, silicone insert across full width of base, additional tray insert, 10" x 15" x 2 5/8" (25 cm x 38 cm x 6.6 cm)	18+
NEW * 2145-03	Silver-gray color, 10.5" x 21" x 2.5" (27 cm x 53 cm x 6.3 cm)	12+ VATS Instruments

All trays are steam, flash, ETO, and plasma sterilization compatible.

Silicone inserts are Latex-Safe.


* Firm Sale

Custom modifications available upon request

LATEX-SAFE



INSTRUMENT CARE PRODUCTS ■ Sterilization Trays ■ Plastic



■ SCANLAN® Plastic Sterilization Trays, Silicone Mat

Cat. No.	Description	Instrument Capacity
2083-01	Metal clasps, 6" x 2" x .7" (15 x 5.1 x 1.8 cm)	1-2
2083-02	Metal clasps, 6" x 3" x .7" (15 x 7.6 x 1.8 cm)	3
2083-03	Metal clasps, 7.8" x 3" x .7" (19.8 x 7.6 x 1.8 cm)	3
2083-05	Metal clasps, 9.6" x 5.5" x .7" (24.4 x 14 x 1.8 cm)	5-10
2081-11	Plastic clasps, 10" x 6" x .7" (25 x 15 x 1.8 cm)	6-10
* 2082-10	Plastic clasps, deep base, 10" x 6" x 1.5" (25 x 15 x 3.8 cm)	6-10
* 2081-21	Plastic clasps, with additional tray insert and extra silicone mat, 10" x 6" x 1.5" (25 x 15 x 3.8 cm)	10+
* 8008-66	Metal clasp, 12" x 9" x .7" (30.5 x 23 cm x 1.8 cm)	8-12
* 8008-62	SCANLAN® Valued Instrument Protection (V.I.P.™) System: tray with metal clasp and soft-sided carrying case, 12" x 9" x .7" (30.5 x 23 cm x 1.8 cm)	8-12

All trays are steam, flash, ETO, and plasma sterilization compatible.






Silicone inserts are Latex-Safe.

* Firm Sale

LATEX-SAFE



INSTRUMENT CARE PRODUCTS ■ Sterilization Trays ■ Plastic



■ SCANLAN® Plastic Sterilization Trays, Silicone Mat

Cat. No.	Description	Instrument Capacity
2081-12	Plastic clasps, 15" x 10" x .7" (38 x 25 x 1.8 cm)	10-15
2082-20	Plastic clasps, deep base, 15" 10" x 1.5" (38 x 25 x 3.8 cm)	10-15
* 2081-23	Plastic clasps, with additional 1/2 tray insert and extra silicone mat, 15" x 10" x 1.5 " (38 x 25 x 3.8 cm)	15+
* 2081-24	Metal clasps, with additional tray insert and extra silicone mat, 15" x 10" x 1.5 " (38 x 25 x 3.8 cm)	20-30
* 2083-16	Metal clasps, 13" x 7.5" x .7 " (33.1 x 19.1 x 1.8 cm)	10-15
* 2084-12	Metal clasps, with additional tray insert and extra silicone mat, 9.5" x 8" x 1.5 " (24.1 x 20.3 x 3.8 cm)	12-15
* 2084-20	Metal clasps, with additional tray insert and extra silicone mat, 17.8" x 9.7" x 1.5 " (45.2 x 25 x 3.8 cm)	24-30+

All trays are steam, flash, ETO, and plasma sterilization compatible.

Silicone inserts are Latex-Safe.

* Firm Sale

LATEX-SAFE

