




**Tubo para Hipospádia**

Tubo para Hipospádia

**Tube for Hypospadias**

Tube pour Hypospadias



**SILIMED** 

**PT**

**Tubo para Hipospádia**

**O Tubo para Hipospádia é um cateter de elastômero de silicone transparente.**

**É indicado na utilização em cirurgias para correção de Hipospádia, visando proporcionar a normalidade estética e funcional do pênis e da uretra. Este produto é indicado para um tratamento com duração aproximada de 15 dias.**

**Fornecido estéril.**

**ES**

**Tubo para Hipospádia**

El Tubo para Hipospadia es un catéter de elastómero de silicona transparente.

Es indicado para uso en cirugías correctivas de Hipospadia, con miras a proporcionar una normalidad estética y funcional al pene y a la uretra. Este producto es indicado para tratamiento de aproximadamente 15 días de duración.

Entregado estéril.

**EN**

**Tube for Hypospadias**

**The Tube for Hypospadias is a transparent silicone elastomer catheter. It is indicated for use in surgeries for correction of Hypospadias, for the purpose of making the penis and urethra aesthetically and functionally normal. This product is indicated for treatment of approximately 15 days.**

**Supplied sterile.**

**FR**

**Tube pour Hypospadias**

Le Tube pour Hypospadias est un Cathéter d'élastomère de silicone transparent.

Il est indiqué en chirurgie de la correction d'Hypospadias et visant à rétablir la normalité esthétique et fonctionnelle du pénis et de l'urètre. Ce produit est indiqué pour un traitement approximatif de 15 jours.

Fourni stérile.

**20**



REF	K >			QTD
	A	ø B		
	mm	French	mm	
3600-006H1	450	6	2,0	1
3600-006H3	450	6	2,0	3
3600-008H1	450	8	2,7	1
3600-008H3	450	8	2,7	3
3600-010H1	450	10	3,3	1
3600-010H3	450	10	3,3	3
3600-012H1	450	12	4,0	1
3600-012H3	450	12	4,0	3
3600-014H1	450	14	4,7	1
3600-014H3	450	14	4,7	3
3600-016H1	450	16	5,3	1
3600-016H3	450	16	5,3	3

